

Anti-Monkey IgG Monoclonal Antibodies

ORDERING INFORMATION

Catalog No.	Clone No.	MAb Subtype	Library Pack No.	100ug/clone
14801	5C12.D4	IgG1	148101	2 clones
14804	4D8.B10	IgG1		

Quantities: 100ug, 500ug, 1mg

Format: Protein G-purified antibody in PBS, pH 7.4.

SPECIFICATION SUMMARY

Antigen: Chromatographically purified Rhesus monkey IgG.

Host Species: Mouse

Specificity: These antibodies recognize Rhesus monkey IgG. They do not cross-react with human, bovine, equine, or ovine IgG.

APPLICATIONS

These antibodies may be used in ELISA to detect and quantitate Rhesus monkey IgG.

DILUTION INSTRUCTIONS

Dilute in PBS or medium that is identical to that used in the assay system.

STORAGE AND STABILITY

These antibodies are stable for at least one (1) year at -20°C to -70°C. Store product in appropriate aliquots to avoid multiple freeze-thaw cycles.

For in vitro investigational use only. Not for use in therapeutic or diagnostic procedures.